

McKesson  
WoundCare Companion™

# Product resource guide

Information about our  
advanced wound care products  
and solutions



# McKesson WoundCare Companion™



**Caring for your patients is a top priority, but advanced wound care can be a challenge.** Let McKesson Medical-Surgical help you customize your treatment plan for wound care, saving staff time while helping your patients heal.

This Product Resource Guide provides details on advanced wound care products and information clinicians need to help them comply with CMS requirements.

With options ranging from trusted national brand manufacturers to our own McKesson Brands portfolio, you can easily find the appropriate products for your patients.

Ask your McKesson Medical-Surgical Account Manager about taking advantage of these opportunities today or visit our full McKesson WoundCare Companion™ program at [mms.mckesson.com/mckesson-woundcare-companion](https://mms.mckesson.com/mckesson-woundcare-companion)



# Table of contents

## Advanced skin care

Swift Medical™ skin and wound solutions	4
Adhesive removers	5
Antifungals	6
Moisture-wicking fabric with antimicrobial silver	6
Moisture barrier ointments	7
Moisturizers	8
No-sting liquid barrier film	9
Skin cleansers	10

## Advanced wound care

Wound measuring guides	11
Wound cleansers	12
Antimicrobial wound cleansers	12

### Wound dressings

Bordered gauze	13
Calcium alginates and fiber gelling dressings	14-16
Collagens with or without antimicrobial silver	17
Composite dressings	18
Silicone contact layer	19
Foam dressings	20-26
Hydrocolloids	26
Hydrogels	27
Hydrophilic wound dressings	27
Medical grade honey	28
Super absorptive dressings	29
Non-adherent dressings	30
Xeroform dressings	30
Silicone tape	31
Surgical tape	32
Transparent films, acrylic and acrylate dressings	33

### Negative pressure wound therapy

Negative pressure wound therapy	34-35
---------------------------------	-------

### Compression dressings

Multi-layer compression bandage systems	36
Tubular support bandages	37
Tubular support retainer dressings	37
Unna boots	38

### Infection prevention

Infection prevention	38
Wound care resources	39



# Swift Medical™ skin and wound solutions

## Setting a new standard

Managing wound care is one of your greatest challenges. McKesson is working with Swift Medical™, a leader in digital wound care management to deliver an accessible, portable and easy to use technology. Swift Skin and Wound is a wound engagement tool that can:

- **Ease your wound care workflows**
- **Standardize high quality of care**
- **Track healing progress for your population**
- **Reduce wound management errors**
- **Reduce the cost of wound care**

**1** With an app that makes capturing wound care information as easy as taking a photo, wound management and tracking is simple and accessible from any mobile device.



**2** Swift offers real-time dashboards to help you identify patients who need immediate attention and make rapid treatment decisions.



## Other benefits include:

- **Increased visibility** into your wound data to empower your team to make better informed and faster care decisions.
- **Real-time collaboration** between clinicians to improve quality and coordination of care.
- **Online and offline capability** to enable clinicians to provide wound care anywhere, regardless of internet connection.
- **Seamless integration with electronic medical records (EMR)** to streamline documentation, reduce errors and improve user satisfaction.

## What has Swift Medical brought to customers?

### Improved quality

**13%** reduction in re-hospitalization rate<sup>1</sup>

Up to **70%** improvement in days to heal a wound<sup>1</sup>

### Increased efficiency

**20%** reduction in visits<sup>1</sup>

**58%** reduction in advanced wound care supply costs<sup>2</sup>

**78%** reduction in wound assessment and documentation time<sup>3</sup>

### Higher staff satisfaction

**85%** overall clinical satisfaction<sup>4</sup>

To find out more, please visit [mms.mckesson.com/swift](https://mms.mckesson.com/swift). To set up a demonstration, please contact your McKesson Account Manager or call **800.654.7240**

1 - [https://2024649.fs1.hubspotusercontent-na1.net/hubfs/2024649/+%20DOCUMENTS%20\(.pdf%20and%20.png\)/+%20Brochures%20and%20One%20Pagers/Swift%20Medical%20Alternate%20Case%20Study.pdf](https://2024649.fs1.hubspotusercontent-na1.net/hubfs/2024649/+%20DOCUMENTS%20(.pdf%20and%20.png)/+%20Brochures%20and%20One%20Pagers/Swift%20Medical%20Alternate%20Case%20Study.pdf)

2 - [https://2024649.fs1.hubspotusercontent-na1.net/hubfs/2024649/+%20DOCUMENTS%20\(.pdf%20and%20.png\)/+%20Brochures%20and%20One%20Pagers/Swift%20Medical%20Cornerstone%20Case%20Study.pdf](https://2024649.fs1.hubspotusercontent-na1.net/hubfs/2024649/+%20DOCUMENTS%20(.pdf%20and%20.png)/+%20Brochures%20and%20One%20Pagers/Swift%20Medical%20Cornerstone%20Case%20Study.pdf)

3 - <https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0271742>

4 - 2021 Swift Medical Customer Survey



# Advanced skin care product portfolio

Maintaining skin integrity is the first step in preventing wounds.

Our products are available in various units of measure – from single products to cases. Contact your McKesson Medical-Surgical Account Manager or check out [mms.mckesson.com](http://mms.mckesson.com) to order today.

## Adhesive removers

Adhesive removers are solvents formulated to dissolve adhesives and assist in removing dressings, skin protectors, tapes and any product that adheres to the skin.



PRODUCT DESCRIPTION

SIZE

MFR #

### 3M™

Remover lotion

4 oz. 20/BX

8610



### MCKESSON

Adhesive remover wipes

4 oz. 50/BX

176-5729



### SMITH & NEPHEW

UNI-SOLVE™ adhesive remover wipes

4 oz. 50/BX

402300

# Antifungals

Antifungals are medicated powders, creams, lotions or ointments used topically to treat fungus and yeast.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MCKESSON</b>		
Thera™ Antifungal Body Powder	3 oz.	53-AFP3
Thera™ Antifungal Body Cream	3 oz.	53-AFC4
<b>SMITH &amp; NEPHEW</b>		
SECURA™ Antifungal Cream	3 ¼ oz.	59432900
SECURA™ Antifungal Greaseless	2 oz.	59432800

# Moisture wicking fabric with antimicrobial silver

Moisture wicking fabric is a soft 100% polyester fabric with ionic silver that wicks and translocates moisture, keeping the skin dry in skin folds and other skin-on-skin contact areas.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>COLOPLAST</b>		
InterDry®	10" x 36" Packet	7912
InterDry®	10" x 144" Roll	7910
<b>MCKESSON</b>		
Moisture Wicking Fabric	10" x 36" Sheet	16-1036S
Moisture Wicking Fabric	10" x 144" Roll	16-10144R

# Moisture barriers

Moisture barriers are creams, ointments or pastes that protect the skin from damage associated with incontinence, perspiration or wound drainage. These products protect the skin, create proper moisture balance and prevent transepidermal water loss.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Cavilon™ Durable Barrier Cream	2 gram Sachet	3353
Cavilon™ Durable Barrier Cream	1 oz.	3354
Cavilon™ Durable Barrier Cream	3 ¼ oz.	3355



<b>MCKESSON</b>		
Thera™ Moisturizing Body Shield	4 gram Packet	53-MS4G
Thera™ Moisturizing Body Shield	4 oz.	53-MS4
Thera™ Dimethicone Body Shield	4 oz.	53-DS4
Thera™ Calazinc Body Shield	4 gram Packet	53-CZ4G
Thera™ Calazinc Body Shield	4 oz.	53-CZ4
Thera™ Silicone Skin Guard	4 oz.	53-SSG4



<b>SMITH &amp; NEPHEW</b>		
SECURA™ Protective Cream	1.75 oz.	59431100
SECURA™ Protective Cream	2.75 oz.	59431200
SECURA™ Extra Protective Cream	3.25 oz.	59432400
SECURA™ Extra Protective Cream	7.7 oz.	59432500
SECURA™ Dimethicone Protectant	4 oz.	59432200
SECURA™ Dimethicone Protectant	3.5 gram Packet	59435000
SECURA™ Protective Ointment	2.47 oz.	59431500
SECURA™ Protective Ointment	5.6 oz.	59431600
SECURA™ Protective Ointment	3.5 gram Packet	59434800



# Moisturizers

Moisturizers help maintain skin integrity by hydrating and soothing dry or irritated skin.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Cavilon™ 3-in-1 Incontinence Care Lotion	8 oz.	3383
Cavilon™ Extra Dry Skin Cream	4 oz.	3386
<b>MCKESSON</b>		
Thera™ Moisturizing Body Cream	1000 mL Refill Bag	53-CRM1000
Thera™ Moisturizing Body Cream	4 oz.	53-CRM4
Thera™ Moisturizing Body Cream	32 oz. Pump Bottle	53-CRM32
<b>SMITH &amp; NEPHEW</b>		
SECURA™ Moisturizing Cream	3 oz.	59431900
SECURA™ Moisturizing Cream	6.5 oz.	5943200

# No-sting liquid barrier film

Liquid barrier films are formulations designed to protect vulnerable areas from the effects of 1) mechanical or chemical injury and 2) excessive moisture due to incontinence, perspiration or wound drainage. They form a transparent film on the skin.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Cavioron™ No Sting Barrier Film	1 mL Wipe	3344
Cavioron™ No Sting Barrier Film	28 mL Bottle	3346
Cavioron™ No Sting Barrier Film	1 mL Wand	3343
Cavioron™ No Sting Barrier Film	3 mL Wand	3345
Cavioron™ Advanced	0.7 mL Wand	5051
Cavioron™ Advanced	2.7 mL Wand	5050



<b>MCKESSON</b>		
No-Sting Skin Barrier Film	Wipe	176-5728



<b>SMITH &amp; NEPHEW</b>		
NO-STING SKIN-PREP™	1 oz. Spray	66800709
NO-STING SKIN-PREP™	3 mL Swab	59420700
NO-STING SKIN-PREP™	1 mL Wipe	59420600

# Rinse-free skin cleansers

Skin cleansers are the first step in the process of cleansing, moisturizing and protecting the skin. Rinse-free formulas clean without water and are pH balanced for optimal skin health.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Cavilon™ No-Rinse Skin Cleanser	8 oz.	3380

## MCKESSON



Rinse-Free Perineal Wash	8 oz.	53-28133
Perineal and Skin Cleanser Rinse-Free	8 oz.	53-28013-8
Foaming Cleanser Rinse-Free	9 oz.	53-22952-9
Thera Foaming Body Cleanser	5 oz. Pump	53-FC5
Thera Foaming Body Cleanser	9 oz. Pump	53-FC9
Thera Moisturizing Body Cleanser	8 oz.	53-MC8
Thera Antimicrobial Body Cleanser	4 oz.	53-AC4
Thera Antimicrobial Body Cleanser	8 oz.	53-A



## SMITH & NEPHEW

SECURA™ Moisturizing Cleanser	4 oz. Spray	59430800
SECURA™ Moisturizing Cleanser	8 oz. Spray	59430900
SECURA™ Total Body Foam Cleanser	4.5 oz.	59430200
SECURA™ Total Body Foam Cleanser	8.5 oz.	59430300




# Advanced wound care product portfolio

Properly selecting and using the appropriate wound care products can be overwhelming. That's why we made it easier to identify effective products and are committed to streamlining wound care.

Our products are available in various units of measure, from single products to cases. Contact your McKesson Medical-Surgical Account Manager or check out [mms.mckesson.com](https://mms.mckesson.com) to order today.

## Wound measuring guides

These measuring guides are disposable and made for single-use.



PRODUCT DESCRIPTION	MFR #
<b>3M™</b> 6" paper ruler	533-LR/50
<b>MCKESSON</b> Plastic circular grid	533-30012100
<b>PURITAN MEDICAL PRODUCTS</b> 6" measuring guide-polystyrene shaft and foam tip	25-1506 1PF DM

# Wound cleansers

Wound cleansers are solutions used to remove debris and exudate from the wound surface.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Wound Cleanser	4 oz.	91101



<b>MCKESSON</b>		
Dermal Wound Cleanser	4 oz.	1718
Dermal Wound Cleanser	8 oz.	1719
Dermal Wound Cleanser	16 oz.	1720
Saline Wound Flush	7 1/10 oz.	37-6507



<b>WINCHESTER LABS</b>		
Saljet® Sterile Saline Solution	30 mL	64938-009-001

# Antimicrobial wound cleansers

These cleansers contain antimicrobial agents if contamination is a consideration.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MCKESSON</b>		
Puracyn® Plus Professional Pump Bottle	4 oz.	186-6505
Puracyn® Plus Professional Pump Bottle	8 1/2 oz.	186-6509
Puracyn® Plus Professional Flip Top Bottle	16 9/10 oz.	186-6516
Puracyn® Plus Professional Spray Bottle	16 9/10 oz.	186-6517
Puracyn® Plus Professional Instill Application	16 9/10 oz.	186-6525
Puracyn® Plus Wound and Skin Care Solution	8 1/2 oz.	186-6002
Puracyn® Plus Wound and Skin Care Solution	16 9/10 oz.	186-6004



<b>SMITH &amp; NEPHEW</b>		
Dermal Wound Cleanser	8 oz.	59449200
Dermal Wound Cleanser	16 oz.	449000

# Bordered gauze

Bordered gauze offers a one-step application and is an alternative to gauze and tape. They are designed with a non-adherent central gauze pad to allow for absorption of light to moderate drainage.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Medipore™ +Pad	1" x 1½"	3562
Medipore™ +Pad	1" x 2 ⅜"	3564
Medipore™ +Pad	1¾" x 2 ⅜"	3566
Medipore™ +Pad	4" x 4"	3568
Medipore™ +Pad	1¾" x 2 ⅜"	3569
Medipore™ +Pad	1¾" x 6"	3570
Medipore™ +Pad	1¾" x 8"	3571
Medipore™ +Pad	1¾" x 11 ¾"	3573



<b>MCKESSON</b>		
Island Dressing	2" x 2"	16-89022
Island Dressing	4" x 4"	16-89044
Island Dressing	4" x 6"	16-89046
Island Dressing	4" x 8"	16-89048
Island Dressing	4" x 10"	16-89041
Island Dressing	4" x 14"	16-89042
Island Dressing	6" x 6"	16-89066
Island Dressing	6" x 8"	16-89068



<b>MOLNLYCKE</b>		
Mepore®	2½" x 3"	670800
Mepore®	3 ⅜" x 4"	670900
Mepore®	3 ⅜" x 6"	671000
Mepore®	3 ⅜" x 8"	671100
Mepore®	3 ⅜" x 10"	671299
Mepore®	3 ⅜" x 12"	671300



<b>SMITH &amp; NEPHEW</b>		
COVRSITE™ Adhesive	4" x 4"	59714100
COVRSITE™ Adhesive	6" x 6"	59714400




# Calcium alginates and other gelling fiber dressings

Calcium alginates and gelling fiber dressings are primary dressings that cover or fill wounds with moderate to heavy drainage.

## Without antimicrobial silver



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Kerracel™ Gelling Fiber	2" x 2"	CWL1032
Kerracel™ Gelling Fiber	4" x 5"	CWL1033
Kerracel™ Gelling Fiber	6" x 6"	CWL1034
Kerracel™ Gelling Fiber	1" x 18" Ribbon	CWL1035



<b>MCKESSON</b>		
Calcium Alginate	2" x 2"	3561
Calcium Alginate	4" x 4 3/4"	3562
Calcium Alginate	4" x 8"	3563
Calcium Alginate	3/4" x 12" Rope	3564
Gelling Fiber	2" x 2"	87200
Gelling Fiber	4" x 4 3/4"	87400
Gelling Fiber	6" x 6"	87600
Gelling Fiber	3/4" x 18" Rope	87180



<b>MOLNLYCKE</b>		
Melgisorb® Plus Calcium Alginate	2" x 2"	252000
Melgisorb® Plus Calcium Alginate	4" x 4"	252200
Melgisorb® Plus Calcium Alginate	4" x 8"	252500
Melgisorb® Plus Calcium Alginate	1 1/4" x 17 3/4" Rope	253500
Exufiber®	2" x 2"	709900
Exufiber®	4" x 5"	709901
Exufiber®	6" x 6"	709903
Exufiber®	0.8" x 17.7" Rope	709909

*Without antimicrobial silver continued on next page*

# Calcium alginates and other gelling fiber dressings *continued*

## Without antimicrobial silver *continued*



PRODUCT DESCRIPTION	SIZE	MFR #
<b>SMITH &amp; NEPHEW</b>		
ALGISITE™ M Calcium Alginate	¾" x 12" Rope	59480400
ALGISITE™ M Calcium Alginate	2" x 2"	59480100
ALGISITE™ M Calcium Alginate	4" x 4"	59480200
ALGISITE™ M Calcium Alginate	6" x 8"	59480300
DURAFIBER Gelling Fiber	2"x 2"	66800559
DURAFIBER Gelling Fiber	4"x 4"	66800560
DURAFIBER Gelling Fiber	4" x 4 ¾"	66800551
DURAFIBER Gelling Fiber	6"x6"	66800561
DURAFIBER Gelling Fiber	¾" x 18" Rope	66800563

## With antimicrobial silver



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Tegaderm™ Alginate Ag	2" x 2"	90301
Tegaderm™ Alginate Ag	4" x 5"	90501
Tegaderm™ Alginate Ag	6" x 6"	90305
Tegaderm™ Alginate Ag	1" x 12" Rope	90307
Kerracel™ Ag Gelling Fiber	2" x 2"	CWL1157
Kerracel™ Ag Gelling Fiber	4" x 5"	CWL1158
Kerracel™ Ag Gelling Fiber	6" x 6"	CWL1159

## MCKESSON



Calcium Alginate with Silver	¾" x 12" Rope	3560
Calcium Alginate with Silver	2" x 2"	3557
Calcium Alginate with Silver	4" x 4 ¾"	3558
Calcium Alginate with Silver	4" x 8"	3559
Gelling Fiber with Silver	2" x 2"	177200
Gelling Fiber with Silver	4" x 4 ¾"	177400
Gelling Fiber with Silver	6" x 6"	177600
Gelling Fiber with Silver	8" x 12"	177812
Gelling Fiber with Silver	¾" x 18" Rope	177180

# Calcium alginates and other gelling fiber dressings *continued*

## With antimicrobial silver



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MOLNLYCKE</b>		
Melgisorb® Ag Calcium Alginate	2" x 2"	255050
Melgisorb® Ag Calcium Alginate	4" x 4"	255100
Melgisorb® Ag Calcium Alginate	6" x 6"	255150
Melgisorb® Ag Calcium Alginate	8" x 12"	255200
Melgisorb® Ag Calcium Alginate	1 2/10" x 18" Rope	255600
Exufiber® Ag Gelling Fiber	2" x 2"	603421
Exufiber® Ag Gelling Fiber	4" x 4 8/10"	603422
Exufiber® Ag Gelling Fiber	4" x 4"	603425
Exufiber® Ag Gelling Fiber	6" x 6"	603423
Exufiber® Ag Gelling Fiber	8" x 12"	603424
Exufiber® Ag Gelling Fiber	4 5/8" x 17 2/3" Rope	603420




<b>SMITH &amp; NEPHEW</b>		
ACTICOAT™ Absorbent Calcium Alginate with Silver	3/4" x 12" Rope	20181
ACTICOAT™ Absorbent Calcium Alginate with Silver	4" x 5"	20381
DURAFIBER Ag Gelling Fiber	2" x 2"	66800570
DURAFIBER Ag Gelling Fiber	4" x 4"	66801174
DURAFIBER Ag Gelling Fiber	4" x 4 3/4"	66800571
DURAFIBER Ag Gelling Fiber	6" x 6"	66800572
DURAFIBER Ag Gelling Fiber	3/4" x 18" Rope	66800574



# Collagens

Collagen wound dressings are derived from a variety of animal sources. They encourage deposition and organization of newly formed collagen fibers and granulation tissue in the wound bed, stimulating new tissue development.



## Without antimicrobial silver

PRODUCT DESCRIPTION SIZE MFR #

### 3M™

PROMOGRAN™ Matrix	4 17/50 sq. in.	PG004
PROMOGRAN™ Matrix	19 1/10 sq. in.	PG019



### MCKESSON

Collagen	2" x 2"	16-1850
Collagen	4" x 4"	16-1852
Collagen Powder	1 gram	1800P



### SMITH & NEPHEW

BIOSTEP™	2" x 2"	66800124
BIOSTEP™	4" x 4"	66800125

## With antimicrobial silver

PRODUCT DESCRIPTION SIZE MFR #

### 3M™

PROMOGRAN™ PRISMA Matrix	4 17/50 sq. in.	MA028
PROMOGRAN™ PRISMA Matrix	19 1/10 sq. in.	MA123
PROMOGRAN™ PRISMA Matrix	19 1/10 sq. in. Rope	MA032



### MCKESSON

Collagen with Silver	2" x 2"	16-1860
Collagen with Silver	4" x 4"	16-1862



### SMITH & NEPHEW

BIOSTEP™ Ag	2" x 2"	66800126
BIOSTEP™ Ag	4" x 4"	66800122



# Composite dressings

Composite dressings are soft, flexible and comfortable for enhanced wearability. They have a water-resistant film layer which acts as an effective barrier to outside contamination and permits showering.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Tegaderm™ +Pad	1" x 1½"	3582
Tegaderm™ +Pad	1" x 2 ⅜"	3584
Tegaderm™ +Pad	4" x 4"	3588
Tegaderm™ +Pad	1 ¾" x 2 ⅜"	3587
Tegaderm™ +Pad	1 ¾" x 4"	3589
Tegaderm™ +Pad	1 ¾" x 6"	3590
Tegaderm™ +Pad	3 ½" x 10"	3591
Tegaderm™ +Pad	3 ½" x 13 ¾"	3593
<b>MCKESSON</b>		
Barrier Island Dressings	4" x 4"	16-89144
Barrier Island Dressings	6" x 6"	16-89166
<b>MOLNLYCKE</b>		
Alldress®	4" x 4"	265329
Alldress®	6" x 6"	265349
Alldress®	6" x 8"	265369
Mepore® Pro	3 ⅜" x 4"	670990
Mepore® Pro	3 ⅜" x 6"	671090
Mepore® Pro	3 ⅜" x 8"	671190
Mepore® Pro	3 ⅜" x 10"	671290
Mepore® Pro	3 ⅜" x 12"	671390
Mepore® Pro	3 ⅜" x 14"	671400
<b>SMITH &amp; NEPHEW</b>		
COVRSITE™ PLUS	4" x 4"	59715000
COVRSITE™ PLUS	6" x 6"	59715100

# Silicone contact layer

Silicone contact layers are dressings with a porous mesh design that allows exudate to pass into a secondary dressing. The gentle silicone adhesive may help reduce pain and trauma upon removal while extending wear time.



## Without antimicrobial silver

PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Adaptic™ Touch Silicone Contact Layer	2" x 3"	500501
Adaptic™ Touch Silicone Contact Layer	3" x 4 1/4"	500502
Adaptic™ Touch Silicone Contact Layer	5" x 6"	500503
Adaptic™ Touch Silicone Contact Layer	8" x 12 3/4"	500504



## MCKESSON

Silicone Contact Layer	3" x 4"	4816
Silicone Contact Layer	4" x 7 2/10"	4817

## MOLNLYCKE

Mepitel® Silicone Contact Layer	2" x 3"	290599
Mepitel® Silicone Contact Layer	3" x 4"	290799
Mepitel® Silicone Contact Layer	4" x 7"	291099
Mepitel® Silicone Contact Layer	8" x 12"	292005
Mepitel® One Silicone Contact Layer	2" x 3"	289100
Mepitel® One Silicone Contact Layer	3" x 4"	289300
Mepitel® One Silicone Contact Layer	4" x 7"	289500
Mepitel® One Silicone Contact Layer	6.8" x 10"	289700
Mepitel® One Silicone Contact Layer	10.8" x 20"	289750



## With antimicrobial silver

PRODUCT DESCRIPTION	SIZE	MFR #
<b>MOLNLYCKE</b>		
Mepitel® Ag Silicone Contact Layer	2" x 3"	39590
Mepitel® Ag Silicone Contact Layer	3" x 4"	390790
Mepitel® Ag Silicone Contact Layer	4" x 5"	390090
Mepitel® Ag Silicone Contact Layer	4" x 7"	391090
Mepitel® Ag Silicone Contact Layer	8" x 12"	392090
Mepitel® Ag Silicone Contact Layer	8" x 20"	392059

# Foam dressings

Foam dressings are used to absorb exudate from a draining wound. They can be bordered or non-bordered, with or without silver, and come in a variety of shapes. Some bordered foams have been shown to reduce friction and shear. Consider these dressings to help protect bony prominences or areas at risk for skin breakdown due to pressure.



## Non-bordered foam

PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Tegaderm™ Foam	2" x 2"	90600
Tegaderm™ Foam	4" x 4"	90601
Tegaderm™ Foam	4" x 8"	90602
Tegaderm™ Foam	8" x 8"	90603
Tegaderm™ Foam	3½" x 3½"	90604
Tegaderm™ Foam, Silicone Face	4" x 4½"	90631
Tegaderm™ Foam, Silicone Face	6" x 6"	90632



<b>MCKESSON</b>		
Foam-Film Backed	4" x 4"	16-4737
Foam-Film Backed	6" x 6"	16-4738
Hydrocellular Foam	4" x 4"	16-4739
Foam, Silicone Face	3" x 3"	4862
Foam, Silicone Face	4" x 4"	4863
Foam, Silicone Face	6" x 6"	4864
Foam, Silicone Face	7" x 7" Sacral	4865
Light Foam, Silicone Face	3" x 3"	4892
Light Foam, Silicone Face	4" x 4"	4893
Light Foam, Silicone Face	6" x 6"	4894

## MOLNLYCKE

Mepilex®, Silicone Face	4" x 4"	294199
Mepilex®, Silicone Face	4" x 5"	294090
Mepilex®, Silicone Face	4" x 8"	294299
Mepilex®, Silicone Face	6" x 6"	294399
Mepilex®, Silicone Face	8" x 8"	294499
Mepilex®, Silicone Face	8" x 20"	294500

*Non-bordered foam continued on next page*

# Foam dressings *continued*

## Non-bordered foam *continued*

PRODUCT DESCRIPTION SIZE MFR #

### MOLNLYCKE *continued*



Mepilex® Heel, Silicone Face	5" x 8"	288100
Mepilex® Lite, Silicone Face	2 2/5" x 3 2/5"	284090
Mepilex® Lite, Silicone Face	4" x 4"	284190
Mepilex® Lite, Silicone Face	6" x 6"	284390
Mepilex® Lite, Silicone Face	8" x 20"	284599
Mepilex® XT, Silicone Face	4" x 4"	211100
Mepilex® XT, Silicone Face	4" x 8"	211200
Mepilex® XT, Silicone Face	6" x 6"	211300
Mepilex® XT, Silicone Face	8" x 8"	211400
Mepilex® Transfer, Silicone Face	6" x 8"	294899
Mepilex® Transfer, Silicone Face	8" x 20"	294599
Lyof foam® Max	3" x 3 2/5"	603200
Lyof foam® Max	4" x 4"	603201
Lyof foam® Max	4" x 8"	603203
Lyof foam® Max	6" x 6"	603204
Lyof foam® Max	6" x 8"	603205
Lyof foam® Max	8" x 8"	603206
Lyof foam® Max T -Fenestrated	3 3/5" x 3 3/5"	603207

### SMITH & NEPHEW



Allevyn Gentle Silicone Face	2" x 2"	66802128
Allevyn Gentle Silicone Face	4" x 4"	66802129
Allevyn Gentle Silicone Face	4" x 8"	66802130
Allevyn Gentle Silicone Face	6" x 6"	66802131
Allevyn Gentle Silicone Face	8" x 8"	66802132
Allevyn Gentle Silicone Face	8" x 20"	66802133
Allevyn Life Silicone Face	2 1/8" x 2 1/8"	66801757
Allevyn Life Silicone Face	4" x 4"	66801748
Allevyn Life Silicone Face	6" x 6"	66801749
Allevyn Life Silicone Face	8 1/4" x 8 1/4"	66801750
Allevyn Life Silicone Face	4" x 8"	66801751
Allevyn Life Silicone Face	8" x 19 3/4"	66801752



# Foam dressings *continued*

## Bordered foam



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Tegaderm™, Acrylate Adhesive	3 ½" x 3 ½"	90610
Tegaderm™, Acrylate Adhesive	4" x 4 ½"	90611
Tegaderm™, Acrylate Adhesive	4" x 4"	90612
Tegaderm™, Acrylate Adhesive	5 ⅝" x 6 ⅙"	90613
Tegaderm™, Acrylate Adhesive	2 ¾" x 3"	90614
Tegaderm™, Acrylate Adhesive	7 ½" x 8 ¾"	90616
Tegaderm™, Acrylate Adhesive	5 ½" x 5 ½"	90619
Tegaderm™, Acrylate Adhesive	2 ¾" x 2 ¾"	90615
Tegaderm™ Foam, Silicone Adhesive	6" x 6"	90646
Tegaderm™ Foam, Silicone Adhesive	6" x 6 ¾"	90647
Tegaderm™ Foam, Silicone Adhesive	7 ¼" x 8 ¾"	90648
Tegaderm™ Foam, Silicone Adhesive	2" x 2"	90643
Tegaderm™ Foam, Silicone Adhesive	3" x 3"	90640
Tegaderm™ Foam, Silicone Adhesive	4" x 4"	90641
Tegaderm™ Foam, Silicone Adhesive	6" x 6"	90642



<b>MCKESSON</b>		
Foam, Acrylic Adhesive	3" x 3"	16-4670
Foam, Acrylic Adhesive	4" x 4"	16-4671
Foam, Acrylic Adhesive	6" x 6"	16-4672
Foam, Acrylic Adhesive	7" x 7"	16-4673
Foam, Silicone Adhesive	3" x 3"	4842
Foam, Silicone Adhesive	4" x 4"	4843
Foam, Silicone Adhesive	6" x 6"	4844
Foam, Silicone Adhesive	7" x 7" Sacral	4845
Light Foam, Silicone Adhesive	3" x 3"	4872
Light Foam, Silicone Adhesive	4" x 4"	4873
Light Foam, Silicone Adhesive	6" x 6"	4874

# Foam dressings *continued*

## Bordered foam

PRODUCT DESCRIPTION

SIZE

MFR #

### MOLNLYCKE



Mepilex® Border Sacrum, Silicone Adhesive	6 <sup>3</sup> / <sub>10</sub> " x 7 <sup>9</sup> / <sub>10</sub> "	282055
Mepilex® Border Sacrum, Silicone Adhesive	8 <sup>7</sup> / <sub>10</sub> " x 9 <sup>8</sup> / <sub>10</sub> "	282455
Mepilex® Border Heel, Silicone Adhesive	8 <sup>7</sup> / <sub>10</sub> " x 9 <sup>1</sup> / <sub>10</sub> "	282790
Mepilex® Border, Silicone Adhesive	3" x 3"	295200
Mepilex® Border, Silicone Adhesive	4" x 4"	295300
Mepilex® Border, Silicone Adhesive	6" x 6"	295400
Mepilex® Border, Silicone Adhesive	6" x 8"	295600
Mepilex® Border, Silicone Adhesive	4" x 8"	295800
Mepilex® Border, Silicone Adhesive	4" x 10"	295850
Mepilex® Border, Silicone Adhesive	4" x 12"	295900
Mepilex® Lite, Silicone Adhesive	1 <sup>6</sup> / <sub>10</sub> " x 2"	581011
Mepilex® Lite, Silicone Adhesive	2" x 5"	581100
Mepilex® Lite, Silicone Adhesive	3" x 3"	581200
Mepilex® Lite, Silicone Adhesive	4" x 4"	581300
Mepilex® Lite, Silicone Adhesive	6" x 6"	581500
Mepilex® Border Flex, Silicone Adhesive	3" x 3"	595200
Mepilex® Border Flex, Silicone Adhesive	4" x 4"	595300
Mepilex® Border Flex, Silicone Adhesive	6" x 6"	595400
Mepilex® Border Flex, Silicone Adhesive	6" x 8"	595600
Mepilex® Border, Silicone Adhesive	4" x 8"	295800
Mepilex® Border, Silicone Adhesive	4" x 10"	295850
Mepilex® Border, Silicone Adhesive	4" x 12"	295900
Mepilex® Post Op, Silicone Adhesive	4" x 6"	496300
Mepilex® Post Op, Silicone Adhesive	4" x 8"	496405
Mepilex® Post Op, Silicone Adhesive	4" x 10"	496455
Mepilex® Post Op, Silicone Adhesive	4" x 12"	496605
Mepilex® Post Op, Silicone Adhesive	4" x 14"	496650

# Foam dressings *continued*

## Bordered foam

PRODUCT DESCRIPTION

SIZE

MFR #



### SMITH & NEPHEW

PRODUCT DESCRIPTION	SIZE	MFR #
Alleyn Life, Silicone Adhesive	4" x 4"	66801067
Alleyn Life, Silicone Adhesive	5 3/4" x 5 3/4"	66801068
Alleyn Life, Silicone Adhesive	6" x 6"	66801069
Alleyn Life, Silicone Adhesive	8 3/4" x 8 3/4"	66801070
Alleyn Life, Silicone Adhesive, Heel	9 4/5" x 9 9/10"	66801304
Alleyn Life, Sacrum, Silicone Adhesive	7" x 7"	66801306
Alleyn Life, Sacrum, Silicone Adhesive	8 1/2" x 9"	66801307
Alleyn Gentle Border Lite, Silicone Adhesive	2" x 2"	66800833
Alleyn Gentle Border Lite, Silicone Adhesive	3" x 3"	66800834
Alleyn Gentle Border Lite, Silicone Adhesive	4" x 4"	66800835
Alleyn Gentle Border Lite, Silicone Adhesive	6" x 6"	66800840
Alleyn Gentle Border	3" x 3"	66800276
Alleyn Gentle Border	4" x 4"	66800270
Alleyn Gentle Border	5" x 5"	66800279
Alleyn Gentle Border	7" x 7"	66800280
Alleyn Gentle Border	6 3/4" x 7 1/16" Multisite	66800959
Alleyn Gentle Border	6 5/8" x 6 3/4" Sacral	66800898
Alleyn Gentle Border	8 1/2" x 9" Sacral	66801031
Alleyn Gentle Border	9" x 9 1/8" Heel	66800506
Alleyn Gentle Border	4" x 8" Surgical	66800900
Alleyn Gentle Border	4" x 10" Surgical	66800264
Alleyn Gentle Border	4" x 12" Surgical	66800265

# Foam dressings *continued*

## Bordered foam, antimicrobial



PRODUCT DESCRIPTION SIZE MFR #

### MCKESSON

Foam with PHMB Acrylic Adhesive	3" x 3"	4973
Foam with PHMB Acrylic Adhesive	4" x 4"	4974
Foam with PHMB Acrylic Adhesive	6" x 6"	4976
Foam with PHMB Acrylic Adhesive	7" x 7" Sacral	4977

### MOLNLYCKE

Mepilex® Border Ag Silicone Adhesive	3" x 3"	395290
Mepilex® Border Ag Silicone Adhesive	4" x 4"	395390
Mepilex® Border Ag Silicone Adhesive	6" x 6"	395490
Mepilex® Border Ag Silicone Adhesive	6" x 8"	395690
Mepilex® Border Ag Silicone Adhesive	4" x 8"	395890
Mepilex® Border Ag Silicone Adhesive	4" x 10"	395790
Mepilex® Border Ag Silicone Adhesive	4" x 12"	395990
Mepilex® Border Post-Op Ag, Silicone Adhesive	4" x 6"	498300
Mepilex® Border Post-Op Ag, Silicone Adhesive	4" x 8"	498400
Mepilex® Border Post-Op Ag, Silicone Adhesive	4" x 10"	498450
Mepilex® Border Post-Op Ag, Silicone Adhesive	4" x 12"	498600
Mepilex® Border Post-Op Ag, Silicone Adhesive	4" x 14"	498650
Mepilex® Border Sacrum Ag, Silicone Adhesive	7 2/10" x 7 2/10"	382090
Mepilex® Border Sacrum Ag, Silicone Adhesive	9 2/10" x 9 2/10"	382490

### SMITH & NEPHEW



Allevyn Ag Gentle, Silicone Adhesive	3" x 3"	66800452
Allevyn Ag Gentle, Silicone Adhesive	5" x 5"	66800453
Allevyn Ag Gentle, Silicone Adhesive	7" x 7"	66800454
Allevyn Ag Sacrum, Silicone Adhesive	6 3/4" x 6 7/8"	66020983
Allevyn Ag Sacrum, Silicone Adhesive	9" x 9"	66020984

# Foam dressings *continued*

## Non-bordered foam, antimicrobial

PRODUCT DESCRIPTION	SIZE	MFR #
<b>MCKESSON</b>		
Foam with PHMB	4" x 4"	16-4940
<b>MOLNLYCKE</b>		
Mepilex® Ag Silicone Face	4" x 4"	287100
Mepilex® Ag Silicone Face	4" x 5"	287090
Mepilex® Ag Silicone Face	4" x 8"	287200
Mepilex® Ag Silicone Face	6" x 6"	287300
Mepilex® Ag Silicone Face	8" x 8"	287400
Mepilex® Ag Silicone Face	8" x 20"	287500
Mepilex® Transfer Ag, Silicone Face	4" x 5"	394190
Mepilex® Transfer Ag, Silicone Face	6" x 8"	394890
Mepilex® Transfer Ag, Silicone Face	8" x 20"	394590



## SMITH & NEPHEW

Allevyn Ag Non-Adhesive	2" x 2"	66020977
Allevyn Ag Non-Adhesive	4" x 4"	66020978
Allevyn Ag Non-Adhesive	6" x 6"	66020980
Allevyn Ag Non-Adhesive	8" x 8"	66020981

# Hydrocolloids

Hydrocolloids are occlusive and water-resistant. They can be used as a primary or secondary dressing, help maintain a moist wound environment and assist in autolytic debridement.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MCKESSON</b>		
Hydrocolloid, Thin	2" x 2"	1882
Hydrocolloid, Thin	4" x 4"	1883
Hydrocolloid, Thin	6" x 6"	1884
Hydrocolloid, Thin Sacral	6" x 7"	1885
Hydrocolloid, Film Backing	2" x 2"	1886
Hydrocolloid, Film Backing	4" x 4"	1887
Hydrocolloid, Film Backing	6" x 7"	1888
Hydrocolloid, Foam Backing	4" x 4"	1889
Hydrocolloid, Foam Backing	6" x 6"	1890

*Hydrocolloids continued on next page*



# Hydrocolloids *continued*



PRODUCT DESCRIPTION	SIZE	MFR #
<b>SMITH &amp; NEPHEW</b>		
REPLICARE™ Thin	2" x 2 3/4"	59484000
REPLICARE™ Thin	3 1/2" x 5 1/2"	59484100
REPLICARE™ Thin	6" x 8"	59484200
REPLICARE™ Ultra	4" x 4"	59484600
REPLICARE™ Ultra	7" x 8"	59484900

# Hydrogels

Used as a primary dressing for wounds with light or no drainage to add or maintain moisture at the wound bed and assist in autolytic debridement. To help protect peri-wound skin from maceration, apply a moisture barrier ointment or barrier film. Available in amorphous, impregnated gauze and sheet forms.

PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Nu-Gel™ Hydrogel	15 gram	MNG415U
<b>MCKESSON</b>		
Hydrogel	1 oz.	1721
Hydrogel	3 oz.	1722
Hydrogel Impregnated Gauze	2" x 2"	61-53022
Hydrogel Impregnated Gauze	4" x 4"	61-53044
Hydrogel Sheet	2" x 2"	61-82022
Hydrogel Sheet	4" x 4"	61-82044
Puracyn® Plus Antimicrobial Hydrogel	3 oz.	186-6542



<b>MOLNLYCKE</b>		
Normlgel® Ag	1 1/2 oz.	350450



<b>SMITH &amp; NEPHEW</b>		
INTRASITE™ Gel	15 gram	66027311
SOLOSITE™ Gel	3 oz.	449600

# Hydrophilic wound dressings

Hydrophilic Wound Dressing with Zinc contains Vitamins A and B6, Calcium and Magnesium. Formulated with a patented zinc-nutrient formulation, and along with its balanced pH technology, provides a moist, mildly acidic environment conducive to wound healing.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MCKESSON</b>		
Hydrophilic Wound Dressing	3 oz.	61-SPD03

# Medical grade honey

MEDIHONEY®'s high osmotic effect pulls fluid from deeper tissues to the wound surface promoting autolytic debridement. The naturally low pH of MEDIHONEY® helps lower the pH of the wound environment, which can impact chronic and stalled wounds.







PRODUCT DESCRIPTION	SIZE	MFR #
<b>MCKESSON</b>		
MEDIHONEY® Gel	½ oz.	31805
MEDIHONEY® Gel	1 ½ oz.	31815
MEDIHONEY® Paste	0.5 oz.	31505
MEDIHONEY® Paste	1.5 oz.	31515
MEDIHONEY® Hydrogel Sheet	2.4" x 2.4"	31622
MEDIHONEY® Hydrogel Sheet	4.3" x 4.3"	31644
MEDIHONEY® Honeycolloid™	2" x 2"	31222
MEDIHONEY® Honeycolloid™	4" x 5"	31245
MEDIHONEY® Calcium Alginate	¾" x 12" Rope	31012
MEDIHONEY® Calcium Alginate	2" x 2"	31022
MEDIHONEY® Calcium Alginate	4" x 5"	31045

MEDIHONEY® is a registered trademark of Medihoney Pty Ltd. and licensed to Derma Sciences, Inc. MEDIHONEY® is not a private label product of McKesson but is available for distribution through McKesson's exclusive brand of Advanced Wound Care products.

# Super absorbent dressings

Super absorbent dressings are highly absorbent and help to decrease the risk of maceration and keep harmful wound components away from the wound and peri-wound skin. They contain a protective outer layer and non-adherent wound contact layer for easy removal.

	PRODUCT DESCRIPTION	SIZE	MFR #
	<b>3M™</b>		
	KerraMax Care® Super-Absorbent Dressing	2" x 2"	PRD500-025
	KerraMax Care® Super-Absorbent Dressing	4" x 4"	PRD500-050
	KerraMax Care® Super-Absorbent Dressing	5" x 6"	PRD500-100
	KerraMax Care® Super-Absorbent Dressing	4" x 9"	PRD500-120
	KerraMax Care® Super-Absorbent Dressing	8" x 9"	PRD500-240
	KerraMax Care® Super-Absorbent Dressing	8" x 12"	PRD500-380
	KerraMax Care® Super-Absorbent Dressing	8" x 20"	PRD500-600
	<b>MCKESSON</b>		
	Super Absorbent Dressing	4" x 5"	61-89545
	Super Absorbent Dressing	6" x 9"	61-89569
	<b>MOLNLYCKE</b>		
	Mextra® Superabsorbent	5" x 5"	610000
	Mextra® Superabsorbent	5" x 7"	610100
	Mextra® Superabsorbent	5" x 9"	610200
	Mextra® Superabsorbent	7" x 9"	610300
	Mextra® Superabsorbent	9" x 11"	610400
	Mextra® Superabsorbent	9" x 13"	610500
	<b>SMITH &amp; NEPHEW</b>		
	EXU-DRY™	4" x 6"	5999004120
	EXU-DRY™	6" x 9"	5999006

# Non-Adherent Dressing

Non-adherent dressings are designed not to stick to the wound bed, resulting in less pain and trauma upon removal.



## Without antimicrobial silver

PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
ADAPTIC™	3" x 3"	2012
ADAPTIC™	3" x 16"	2014
ADAPTIC™	3" x 8"	2015
ADAPTIC™	5" x 9"	2016

## MCKESSON

Oil Emulsion	3" x 3"	61-77041
Oil Emulsion	3" x 8"	61-77042

## With antimicrobial silver

PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
SILVERCEL™ NON-ADHERENT Antimicrobial Dressing	2" x 2"	900202
SILVERCEL™ NON-ADHERENT Antimicrobial Dressing	4 ¼" x 4 ¼"	900404
SILVERCEL™ NON-ADHERENT Antimicrobial Dressing	4" x 8"	900408
SILVERCEL™ NON-ADHERENT Antimicrobial Dressing	1" x 12"	900112

# Xeroform dressing


Primary dressing that contains 3% Bismuth Tribromophenate to provide bacteriostatic properties and help reduce risk of infection and deodorize the wound. Helps maintain a moist wound healing environment to assist in autolytic debridement for wounds with light or no drainage.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MCKESSON</b>		
Xeroform with Petrolatum	1" x 8"	2201
Xeroform with Petrolatum	2" x 2"	2202
Xeroform with Petrolatum	4" x 4"	2206
Xeroform with Petrolatum	5" x 9"	2207

# Silicone tape

Silicone Tape is ideal for minimizing pain and trauma for sensitive or at-risk skin. May be appropriate for individuals with adhesive allergies.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b> Micropore™ S Silicone Tape	1" x 5 ½ yd	2770-1
Micropore™ S Silicone Tape	2" x 1 ½ yd	2770S-2
Micropore™ S Silicone Tape	2" x 5 ½ yd	2770-2
Micropore™ S Silicone Tape	1" x 1 ½ yd	2770S-1



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MCKESSON</b> Silicone Adhesive Surgical Tape	1" x 5 ½ yd	16-48410
Silicone Adhesive Surgical Tape	2" x 5 ½ yd	16-48420



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MOLNLYCKE</b> Mepitac® Soft Silicone Tape	¾" x 118"	298300
Mepitac® Soft Silicone Tape	1 ½" x 59"	298400



# Surgical tape

Tapes use adhesive to secure wound dressings or other devices.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Medipore™ H Cloth Surgical Tape	1" x 10 yd	2861
Medipore™ H Cloth Surgical Tape	2" x 10 yd	2862
Medipore™ H Cloth Surgical Tape	3" x 10 yd	2863
Medipore™ H Cloth Surgical Tape	4" x 10 yd	2864
Medipore™ H Cloth Surgical Tape	6" x 10 yd	2866
Medipore™ H Cloth Surgical Tape	8" x 10 yd	2868
Medipore™ H Cloth Surgical Tape, single-patient use roll	1" x 2 yd	2860S-1
Medipore™ H Cloth Surgical Tape, single-patient use roll	2" x 2 yd	2860S-2
Medipore™ H Cloth Surgical Tape, single-patient use roll	4" x 2 yd	2860S-4
Medipore™ H Cloth Surgical Tape, single-patient use roll	6" x 2 yd	2860S-6
Durapore™ Surgical Tape	½" x 10 yd	1538-0
Durapore™ Surgical Tape	1" x 10 yd	1538-1
Durapore™ Surgical Tape	2" x 10 yd	1538-2
Durapore™ Surgical Tape	3" x 10 yd	1538-3
Durapore™ Surgical Tape, single-patient use roll	1 1/10" x 1 ½ yd	1538S-1
Durapore™ Surgical Tape, single-patient use roll	2" x 1 ½ yd	1538S-2



<b>MCKESSON</b>		
Premium Transparent Surgical Tape	1" x 10 yd	100197
Premium Transparent Surgical Tape	2" x 10 yd	100198
Transparent Surgical Tape	1" x 10 yd	16-47210
Transparent Surgical Tape	2" x 10 yd	16-47220

# Transparent films, acrylic and acrylate dressings

Transparent films and clear acrylic dressings are used as primary or secondary dressings. They help maintain a moist wound environment. These dressings are water-resistant and transparent.



## PRODUCT DESCRIPTION

## SIZE

## MFR #

### 3M™

Tegaderm™ Transparent Film	2 3/8" x 2 3/4"	1624W
Tegaderm™ Transparent Film	4" x 4 3/4"	1626W
Tegaderm™ Absorbent	4 3/8" x 5"	90801
Tegaderm™ Absorbent	3" x 3 3/4"	90800
Tegaderm™ Absorbent	5 5/8" x 6 1/4"	90803
Tegaderm™ Absorbent	6 5/8" x 7 1/2"	90807
Tegaderm™ Absorbent	5 7/8" x 6"	90802
Tegaderm™ Absorbent	7 7/8" x 8"	90805
Tegaderm™ Absorbent	3 1/2" x 10"	90815



### MCKESSON

Transparent Film, Frame Delivery	2 3/8" x 2 3/4"	4985
Transparent Film, Frame Delivery	4" x 4 3/4"	4986
Transparent Film, Frame Delivery	6" x 8"	4987



### MOLNLYCKE

Mepore® Film	2 2/5" x 2 3/5"	270600
Mepore® Film	4" x 5"	271500
Mepore® Film	4" x 10"	272500
Mepore® Film	6" x 8 1/2"	273000




### SMITH & NEPHEW

OPSITE FLEXIGRID™	2 3/8" x 2 3/4"	66024628
OPSITE FLEXIGRID™	4" x 4 3/4"	66024630
OPSITE FLEXIGRID™	6" x 8"	66024631

# Negative pressure wound therapy

Negative pressure wound therapy (NPWT) involves the application of controlled levels of sub-atmospheric (negative) pressure to a wound. NPWT may stimulate formation of granulation tissue, increases blood flow, may promote angiogenesis, reduces edema, improves wound contraction and protects from external contamination.<sup>1</sup> All devices listed below are considered disposable.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
SNAP™ Dressing Kit	10 cm x 10 cm	SKTF10X10
SNAP™ Dressing Kit	15 cm x 15 cm	SKTF15X15
SNAP™ Dressing Kit	20 cm x 20 cm	SKTF20X20
SNAP™ Therapy Cartridge	60 cc-75 mmHg	SNPA075US
SNAP™ Therapy Cartridge	60 cc-100 mmHg	SNPA100US
SNAP™ Therapy Cartridge	60 cc-125 mmHg	SNPA125US
SNAP PLUS™ Therapy Cartridge	150 cc-125 mmHg	SNPA125PLUS
PREVENA™ PEEL & PLACE™ Negative Pressure Wound Therapy Kit	13 cm	PRE1101US
PREVENA™ PEEL & PLACE™ Negative Pressure Wound Therapy Kit	20 cm	PRE1001US
PREVENA PLUS™ PEEL & PLACE™ Negative Pressure Wound Therapy Kit	35 cm	PRE3201US
PREVENA™ Canister (5 per case)	45 mL	PRE1095
PREVENA PLUS™ Canister	150 mL	PRE4095

1 - Nie B, Yue B. Biological effects and clinical application of negative pressure wound therapy: a review. J Wound Care. 2016 Nov 2;25(11):617-626. doi: 10.12968/jowc.2016.25.11.617. PMID: 27827278.

# Negative pressure wound therapy *continued*



PRODUCT DESCRIPTION	SIZE	MFR #
<b>SMITH &amp; NEPHEW</b>		
PICO 7 - Multi-site-1 Dressing Kit	15 x 20 cm	66022010
PICO 7 - Multi-site-1 Dressing Kit	20 x 20 cm	66022011
PICO 7 - 1 Dressing Kit	10 x 20 cm	66022012
PICO 7 - 1 Dressing Kit	10 x 30 cm	66022013
PICO 7 - 1 Dressing Kit	10 x 40 cm	66022014
PICO 7 - 1 Dressing Kit	15 x 15 cm	66022015
PICO 7 - 1 Dressing Kit	15 x 20 cm	66022016
PICO 7 - 1 Dressing Kit	15 x 30 cm	66022017
PICO 7 - 1 Dressing Kit	20 x 20 cm	66022018
PICO 7 - 1 Dressing Kit	25 x 25 cm	66022019
PICO 7 Multi-site Small - 2 Dressing Kit	15 x 20 cm	66022000
PICO 7 Multi-site Large - 2 Dressing Kit	15 x 30 cm	66022001
PICO 7 - 2 Dressing Kit	10 x 20 cm	66022002
PICO 7 - 2 Dressing Kit	10 x 30 cm	66022003
PICO 7 - 2 Dressing Kit	10 x 40 cm	66022004
PICO 7 - 2 Dressing Kit	15 x 15 cm	66022005
PICO 7 - 2 Dressing Kit	15 x 20 cm	66022006
PICO 7 - 2 Dressing Kit	15 x 30 cm	66022007
PICO 7 - 2 Dressing Kit	20 x 20 cm	66022008
PICO 7 - 2 Dressing Kit	25 x 25 cm	66022009
PICO 7Y, Wound Kit	N/A	66022031
PICO 14 - 2 Dressing Kit	10 x 20 cm	66022042
PICO 14 - 2 Dressing Kit	10 x 30 cm	66022043
PICO 14 - 2 Dressing Kit	10 x 40 cm	66022044
PICO 14 - 2 Dressing Kit	15 x 15 cm	66022045
PICO 14 - 2 Dressing Kit	15 x 20 cm	66022046
PICO 14 - 2 Dressing Kit	15 x 30 cm	66022047
PICO 14 - 2 Dressing Kit	25 x 25 cm	66022049

# Multi-layer compression bandage systems

Multi-layer compression bandage systems are used in venous disease and edema management. These systems provide graduated, sustained compression.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Coban™ 2 Two-Layer Compression System with Stocking	2 <sup>9</sup> / <sub>10</sub> yd X 4" / 4" X 5 <sup>1</sup> / <sub>10</sub> yd	2094N
Coban™ 2 Two-Layer Compression System	6" x 3 <sup>4</sup> / <sub>5</sub> yd / 6" x 4 <sup>9</sup> / <sub>10</sub> yd	20096
Coban™ 2 Two-Layer Compression System	2" x 1 <sup>3</sup> / <sub>10</sub> yd / 2" x 3 yd	2092
Coban™ 2 Two-Layer Compression System	4" x 3 <sup>4</sup> / <sub>5</sub> yd / 4" x 6 <sup>3</sup> / <sub>10</sub> yd	2094XL
Coban™ 2 Lite Two-Layer Compression System with Stocking	4" x 2 <sup>9</sup> / <sub>10</sub> yd / 4" x 5 <sup>1</sup> / <sub>10</sub> yd	2794N



<b>MCKESSON</b>		
Compression Bandage System, 2 Layer	N/A	2006
Compression Bandage System, 4 Layer	N/A	1008

<b>SMITH &amp; NEPHEW</b>		
PROFORE™ Lite Compression Bandaging System, 3 layer	N/A	66000771
PROFORE™ Compression Bandaging System, 4 layer	N/A	66020016



# Tubular support bandage

These are tubular elastic bandages that provide even compression and support.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MCKESSON</b>		
Elastic Tubular Bandage	Med Arm/Sm Ankle	182-13112C
Elastic Tubular Bandage	Lg Arm/Med Ankle/Sm Knee	182-13113D
Elastic Tubular Bandage	Lg Ankle/Med Knee/Sm Thigh	182-13114E
Elastic Tubular Bandage	Lg Knee/Med Thigh	182-13115F
Elastic Tubular Bandage	Lg Thigh	182-13116G
Elastic Tubular Bandage	Sm Trunks	182-13117J



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MOLNLYCKE</b>		
Tubigrip®	Med Arm/Sm Ankle	1443
Tubigrip®	Lg Arm/Med Ankle/Sm Knee	1437
Tubigrip®	Knee/Sm Thigh	1434
Tubigrip®	Lg Knee/Med Thigh	1438
Tubigrip®	Lg Thigh	1439
Tubigrip®	Sm Trunks	1440

# Tubular net retainer dressing

Tubular net retainer dressings are secondary dressings which can keep bandages and devices securely in place without adhesives.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>3M™</b>		
Tubular Elastic Retainer Net	Size 1	MSVP114701
Tubular Elastic Retainer Net	Size 2	MSVP114702
Tubular Elastic Retainer Net	Size 3	MSVP114703
Tubular Elastic Retainer Net	Size 4	MSVP114704
Tubular Elastic Retainer Net	Size 5	MSVP114705
Tubular Elastic Retainer Net	Size 6	MSVP114706
Tubular Elastic Retainer Net	Size 7	MSVP114707
Tubular Elastic Retainer Net	Size 8	MSVP114708
Tubular Elastic Retainer Net	Size 9	MSVP114709
Tubular Elastic Retainer Net	Size 10	MSVP114710
Tubular Elastic Retainer Net	Size 11	MSVP114711

# Unna boots

Unna boots are non-hardening zinc impregnated bandages used for compression therapy. They are available with or without calamine.



PRODUCT DESCRIPTION	SIZE	MFR #
<b>MCKESSON</b>		
Unna Boot with Zinc Oxide, Gelatin Free	3" x 10 yd	2066
Unna Boot with Zinc Oxide, Gelatin Free	4" x 10 yd	2067
Unna Boot with Zinc Oxide and Calamine, Gelatin Free	3" x 10 yd	2068
Unna Boot with Zinc Oxide and Calamine, Gelatin Free	4" x 10 yd	2069
<b>SMITH &amp; NEPHEW</b>		
VISCOPASTE™ PB7 Zinc Paste Bandage	3"	4956

# Infection prevention

Infection prevention protocols help to protect the patient and the wound care provider. The Centers for Medicare and Medicaid Services (CMS) requires that long term care facilities and home health agencies maintain an infection prevention program based upon CDC or other nationally recognized infection control guidelines. Make the commitment to help prevent infections by utilizing safe and healthy practices every day. You have the power to be an infection preventionist and we are here to help.

Check out our infection prevention program at: [mms.mckesson.com/uprevent](http://mms.mckesson.com/uprevent)



# Wound care resources

In addition to the McKesson WoundCare Companion™ program, the resources below are available for additional information on wound care education, prevention and treatment.



**3M**  
[www.3m.com](http://www.3m.com)

---



**Association for the Advancement of Wound Care**  
[aawconline.org](http://aawconline.org)

---



**Center for Medicare and Medicaid Services**  
[www.CMS.gov](http://www.CMS.gov)

---



**Coloplast**  
[www.coloplast.us](http://www.coloplast.us)

---



**Innovacyn**  
[www.innovacyn.com](http://www.innovacyn.com)

---



**Kestrel Wound Source**  
[www.woundsource.com](http://www.woundsource.com)

---



**McKesson Brands**  
[mms.mckesson.com/mckesson-brand](http://mms.mckesson.com/mckesson-brand)

---



**McKesson Training and Compliance**  
[/mms.mckesson.com/content/our-services-solutions/training-compliance/](http://mms.mckesson.com/content/our-services-solutions/training-compliance/)

---



**Mölnlycke**  
[www.molnlycke.us](http://www.molnlycke.us)

---



**National Pressure Injury Advisory Panel**  
[www.npiap.org](http://www.npiap.org)

---



**Smith & Nephew**  
[www.smith-nephew.com](http://www.smith-nephew.com)

---



**Wound Care Education Institute**  
[www.wcei.net](http://www.wcei.net)

---



**Wound, Ostomy, Contenance Nurses Society**  
[www.wocn.org](http://www.wocn.org)

---

## Have clinical questions?

McKesson Medical-Surgical's certified wound specialists are ready to help you today.

**Call our Clinicians at 877.611.0081**

## Need educational support?

McKesson Academy™ delivers clinical training and education options so your staff can take CE courses anytime, anywhere.

**For more information, contact your Account Manager or call our Clinical Resource Team at 877.611.0081**

## To learn more about WoundCare Companion™

**visit [mms.mckesson.com/mckesson-woundcare-companion](https://mms.mckesson.com/mckesson-woundcare-companion)**

or call your McKesson Medical-Surgical Account Manager.

Be advised that information contained herein is intended to serve as a useful reference for informational purposes only and is not complete clinical information. This information is intended for use only by competent healthcare professionals exercising judgment in providing care. McKesson Medical-Surgical cannot be held responsible for the continued currency of or for any errors or omissions in the information.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers, not McKesson Medical-Surgical or its affiliates ("McKesson"), and have not been independently verified by McKesson. McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use. Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson Medical-Surgical is providing HCPCS codes based upon information provided by the product manufacturers. McKesson Medical-Surgical does not make representations about the accuracy of this information. The HCPCS codes do not indicate coverage or payment under any federal or state healthcare program. The entity that submits claims to any federal or state healthcare program has the responsibility to verify the accuracy of HCPCS codes it submits for payment purposes.

All trademarks and registered trademarks are the property of their respective owners.

© 2023 McKesson Medical-Surgical Inc. 2022\_1700562