



- Designed to make access clear and easy.
- Available in convenient lengths and configurations.
- Ideal for Surgical Centers, Oncology Centers, Emergency Medicine Services and Physician Offices where 24-hour use is applicable.

SafeDAY™ IV Administration and Extension Sets

Affordable, IV Sets with Needleless Connectors
limited to 24-hour maximum use

Easy to:

Visually Inspect



Clear design helps verify proper flush.

Clean



Smooth surfaces permit easy cleaning.

Grip



Easy grip barrel designed to minimize slips and touch contamination.







Connect



Open fluid path designed to minimize areas where bacteria can collect.











Ordering Information

SafeDAY™ sets are restricted for 24 hour maximum use.

| REF No. | Description | Units/Case |
|---------|---|------------|
| 415104 |  SafeDAY (LAD) valve for aspiration, injection or gravity flow of fluid upon insertion of a male luer fitting. Not made with DEHP or natural rubber latex. Priming volume: 0.22 mL Power injection rated to 400 PSI and 15mL/sec. | 100 |
| 352290 |  Standard bore extension set with bonded SafeDAY (LAD) valve, non-removable slide clamp, SPIN-LOCK® connector. Not made with natural rubber latex or DEHP. Priming volume: 1.6 mL, Length: 8 in. (20 cm) | 100 |
| 352905 |  Small bore extension set with bonded SafeDAY (LAD) valve, non-removable purple slide clamp, SPIN-LOCK connector. Not made with natural rubber latex or DEHP. Priming volume: 0.6 mL, Length: 8 in. (20 cm). Power injection rated to 300 psi and 10 mL/sec. | 100 |
| 352906 |  Small bore extension set with bonded SafeDAY (LAD) valve, non-removable slide clamp, SPIN-LOCK connector. Not made with natural rubber latex or DEHP. Priming volume: 0.7 mL, Length: 14 in. (36 cm) | 100 |
| 352633 |  Standard bore extension set with SafeLine® injection site 4 in. above distal end, non-removable slide clamp, SafeDAY valve 2 in. above distal end, sliding SPIN-LOCK connector. Not made with natural rubber latex or DEHP. Priming volume: 1.6 mL, Length: 9 in. (23 cm) | 50 |
| 352644 |  Standard bore extension set with female luer lock connector, SafeDAY valves 25 in. and 6 in. above distal end, sliding SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 4.9 mL, Length: 35 in. (89 cm) | 50 |

Ordering Information





SafeDAY™ sets are limited to 24-hour maximum use.

| REF No. | Description | Units/Case |
|---------|---|------------|
| 352641 |  (15 drops/mL) Universal spike, roller clamp, slide clamp, SafeDAY valve 8 in. above distal end, sliding SPIN-LOCK® connector. Not made with DEHP or natural rubber latex. Priming volume: 17 mL, Length: 86 in. (218 cm) | 50 |
| 352636 |  (15 drops/mL) Universal spike, slide clamp, backcheck valve, SafeDAY valves 87 in. and 14 in. above distal end, roller clamp, sliding SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 18 mL, Length: 105 in. (267 cm) | 50 |
| 352639 |  (15 drops/mL) Universal spike, backcheck valve, SafeDAY valves 69 in. and 6 in. above distal end, roller clamp, sliding SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 15 mL, Length: 84 in. (213 cm) | 50 |
| 352645 |  (15 drops/mL) Universal spike, roller clamp, SafeDAY valves 10 in. and 4 in. above distal end, sliding SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 17 mL, Length: 88 in. (224 cm) | 50 |
| 352640 |  (15 drops/mL) Universal spike, backcheck valve, SafeDAY valves 85 in., 41 in. and 6 in. above distal end, roller clamp, slide clamp, sliding SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 20 mL, Length: 112 in. (284 cm) | 50 |
| 352646 |  (15 drops/mL) Universal spike, backcheck valve, SafeDAY valves 100 in. and 8 in. above distal end, removable high-flow stopcock, sliding SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 23 mL, Length: 125 in. (318 cm) | 50 |
| 352643 |  (15 drops/mL) Universal spike, backcheck valve, SafeDAY valves 76 in. and 6 in. above distal end, roller clamp, slide clamp, Supor® Membrane 0.2 micron air-eliminating filter, sliding SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 22 mL, Length: 104 in. (264 cm) | 50 |
| 362420 |  Infusomat® Space® Pump Set, (15 drops/mL) Universal spike, slide clamp, pressure limited check valve, SafeDAY injection sites 85 in. and 6 in. above distal end, free flow protection clamp, backcheck valve at distal Y-site, roller clamp, slide clamp, SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 22 mL, Length: 120 in. (305 cm) | 24 |
| 352642 |  (15 drops/mL) Universal spike, roller clamp, SafeDAY valves 20 in. and 14 in. above distal end, and removable 8 in. extension set with slide clamp, sliding SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 19 mL, Length: 97 in. (246 cm) | 50 |
| 352635 |  (60 drops/mL) Universal spike, slide clamp, backcheck valve, SafeDAY valves 86 in. and 14 in. above distal end, sliding SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 18 mL, Length: 105 in. (267 cm) | 50 |

To order, call your B. Braun Representative at 800-227-2862.

Ordering Information

SafeDAY™ sets are limited to 24-hour maximum use.

| REF No. | Description | Units/ Case |
|---------|---|----------------|
| 352630 |  <p>(10 drops/mL) Universal spike with SafeLine® injection site 12 in. above distal end, SafeDAY valve 6 in. above distal end, sliding SPIN-LOCK connector. Not made with natural rubber latex or DEHP. Priming Volume: 16 mL, Length: 88 in. (224cm)</p> | 50 |
| 352631 |  <p>(10 drops/mL) Universal spike, roller clamp, SafeLine injection site 16 in. above distal end, SafeDAY valve 10 in. above distal end, and removable extension set with slide clamp, sliding SPIN-LOCK connector. Not made with natural rubber latex or DEHP. Priming Volume: 15 mL, Length: 88 in. (224cm)</p> | 50 |
| 352632 |  <p>(10 drops/mL) Universal spike, roller clamp, SafeDAY valve 16 in. above distal end, SafeLine injection site 10 in. above distal end, and removable extension set with SafeDAY valve 4 in. above distal end, sliding SPIN-LOCK connector. Not made with natural rubber latex or DEHP. Priming Volume: 15 mL, Length: 84 in. (213cm)</p> | 50 |
| 352634 |  <p>(10 drops/mL) Universal spike, roller clamp, two SafeDAY valves 12 in. and 6 in. above distal end, sliding SPIN-LOCK connector. Not made with natural rubber latex or DEHP. Priming Volume: 16 mL, Length: 88 in. (224 cm)</p> | 50 |

B. Braun Medical Inc. | Bethlehem PA | 800-227-2862 | bbraunusa.com